

Dimethyl-Histone H3-K9 Rabbit Polyclonal Antibody (Q16695/P68431) For Epigenetic Research | Wb, Chip, If/Icc

Item Number: CM0000390



Introduction

DiMethyl-Histone H3-K9 rabbit polyclonal antibody (Q16695/P68431) validated for WB, ChIP, IF/ICC, and more. Reacts with human, mouse, rat, and other species. Ideal for epigenetic and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H3-K9 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K9
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag