

Dimethyl-Histone H4-K20 (P62805) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000292



Introduction

DiMethyl-Histone H4-K20 (P62805) rabbit monoclonal antibody validated for WB, DB, IHC-P, IF/ICC, and ELISA. Targets dimethylated lysine 20 of histone H4. Cross-reacts with human, mouse, and rat. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H4-K20 Rabbit mAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; DiMethyl-Histone H4-K20
Species	Human
GeneID (human)	8359/8370
GeneID	8359
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11 kDa
Shipping	Ice bag