


Dimethyl-Histone H4-K20 Rabbit Monoclonal Antibody (P62805) - Validated For Wb, Db, Ihc-P, Elisa In Epigenetics Studies

Item Number: CM0000321



Introduction

Specifically detects dimethylated lysine 20 of human Histone H4 (P62805). Ideal for WB, DB, IHC-P, ELISA with cross-reactivity to Human, Mouse, Rat and predicted wide range. Suitable for chromatin biology and epigenetic IVD applications.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H4-K20 Rabbit mAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; DiMethyl-Histone H4-K20
Species	Human
GeneID (human)	8359/8370
GeneID	8359
Immunogen	A synthetic dimethylated peptide around K20 of human Histone H4 (NP_003539.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag