


Dimethyl-Histone H4-K20 Rabbit Pab (P62805) For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000381



Introduction

Rabbit polyclonal antibody targeting dimethylated histone H4 at lysine 20 (H4K20me2). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat, and a wide range of species. Detects an ~11 kDa band; supplied liquid, shipped on ice bags. Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	DiMethyl-Histone H4-K20 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; DiMethyl-Histone H4-K20
Species	Human
GeneID (human)	H4
GeneID	8359
Immunogen	Synthetic peptide: A synthetic dimethylated peptide around K20 of human histone H4 (NP_003529.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag