

Dlat Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q8Bmf4

Item Number: CM0025726

Introduction

DLAT rabbit polyclonal antibody against mouse dihydrolipoamide acetyltransferase (PDC-E2). Validated in WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. For metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	DLAT Rabbit pAb
Remarks/Alias	DLTA; PDC-E2; 6332404G05Rik
Species	Mouse
GeneID (Human)	1737
GeneID	235339
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 86-305 of mouse DLAT (NP_663589.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8BMF4
Protein Weight	68 kDa
Shipping	Ice bag