

DLx5 Rabbit Mab For Wb, Elisa - P56178

Item Number: CM0009354



Introduction

DLX5 rabbit monoclonal antibody detects human DLX5, a transcriptional factor essential for osteoblast differentiation and bone development. Applications: WB, ELISA. Ideal for bone biology and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	DLX5 Rabbit mAb
Remarks/Alias	SHFM1; SHFM1D; DLX5
Species	Human
Gene ID (Human)	1749
Gene ID	1749
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 150-250 of human DLX5 (P56178).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56178
Protein Weight	32kDa
Shipping	Ice bag