

Dmgdh Rabbit Mab For Wb, If-P, Elisa - Q9Ui17

Item Number: CM0007563



Introduction

DMGDH Rabbit mAb is a monoclonal antibody validated for WB, IF-P, ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying mitochondrial dimethylglycine dehydrogenase function in glycine metabolism.

[Learn More](#)

Parameter	Value
Product Name	DMGDH Rabbit mAb
Remarks/Alias	DMGDHD; ME2GLYDH
Species	Human
Gene ID (Human)	29958
Gene ID	29958
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 472-866 of human DMGDH (NP_037523.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Ui17
Protein Weight	45KD/97kDa
Shipping	Ice bag