

Dnaja4 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q8Ww22

Item Number: CM0003573

Introduction

DNAJA4 Rabbit Monoclonal Antibody validated for WB, IP, ELISA with cross-reactivity to mouse and rat. Ideal for research on membrane-localized chaperone activity, alternative splicing, and zinc-finger protein studies.

[Learn More](#)

Parameter	Value
Product Name	DNAJA4 Rabbit mAb
Remarks/Alias	MST104, MSTP104, PRO1472, DNAJA4
Species	Human
Gene ID (Human)	55466
Gene ID	55466
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-397 of human DNAJA4 (NP_057538.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8WW22
Protein Weight	45kDa
Shipping	Ice bag