

Dnajc17 Rabbit Pab - Q9Nvm6

Item Number: CM0024237



Introduction

Rabbit polyclonal antibody targeting human DNAJC17, suitable for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Recognizes a ~35 kDa protein. Useful for investigating PAX8-mediated transcriptional regulation and chaperone-related research.

[Learn More](#)

Parameter	Value
Product Name	DNAJC17 Rabbit pAb
Remarks/Alias	DNAJC17
Species	Human
GeneID (Human)	55192
GeneID	55192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human DNAJC17 (NP_060633.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVM6
Protein Weight	35kDa
Shipping	Ice bag