

Dnmt1 Rabbit Pab For Wb, If/Icc, Ip, Elisa - P26358

Item Number: CM0028888



Introduction

DNMT1 Rabbit pAb is a rabbit polyclonal antibody against human DNMT1, validated for WB, IF/ICC, IP, ELISA. Ideal for DNA methylation, epigenetic inheritance, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	DNMT1 Rabbit pAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsal; DNMT1
Species	Human
Gene ID (Human)	1786
Gene ID	1786
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 760-1109 of DNMT1 (NP_001370.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag