


Dnmt3A Rabbit Pab For Wb, Ip, Elisa - Q9Y6K1

Item Number: CM0021398



Introduction

High-affinity rabbit polyclonal antibody targeting DNMT3A, validated for Western blot, immunoprecipitation, and ELISA. Detects human, mouse, and rat DNMT3A. Suitable for de novo DNA methylation and epigenetic studies.

[Learn More](#)

Parameter	Value
Product Name	DNMT3A Rabbit pAb
Remarks / Alias	TBR5; HESJAS; DNMT3A2; M.HsaIIIA; DNMT3A
Species	Human
Gene ID (Human)	1788
Gene ID	1788
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-132 of human DNMT3A (NP_072046.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6K1
Protein Weight	102kDa
Shipping	Ice bag