


## Dopamine Transporter/Dat Rabbit Mab For Wb, Ihc-P, Elisa - Q01959

Item Number: CM0005396



### Introduction

Rabbit monoclonal antibody targeting human Dopamine Transporter (DAT) for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying dopaminergic signaling and neurodegeneration. SWISS: Q01959.

[Learn More](#)

| Parameter        | Value                               |
|------------------|-------------------------------------|
| Product Name     | Dopamine Transporter/DAT Rabbit mAb |
| Remarks/Alias    | DAT; DAT1; PKDYS; PKDYS1            |
| Species          | Human                               |
| GeneID (Human)   | 6531                                |
| GeneID           | 6531                                |
| Immunogen        | Synthetic peptide                   |
| Source           | Rabbit                              |
| Category         | Monoclonal Antibodies               |
| Application      | WB, IHC-P, ELISA                    |
| Cross Reactivity | Mouse, Rat                          |
| SWISS            | Q01959                              |
| Protein Weight   | 68kDa                               |
| Shipping         | Ice bag                             |