

Dpp3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ny33

Item Number: CM0030392



Introduction

DPP3 rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, and rat DPP3, a cytosolic zinc metalloprotease that cleaves bioactive peptides like angiotensins and enkephalins. Ideal for IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	DPP3 Rabbit pAb
Remarks/Alias	DPPIII; DPP3
Species	Human
Gene ID	10072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human DPP3 (NP_569710.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NY33
Protein Weight	83kDa
Shipping	Ice bag