


Drd2 Rabbit Polyclonal Antibody For Wb And Elisa - P14416

Item Number: CM0016640



Introduction

Rabbit polyclonal antibody targeting human dopamine D2 receptor (DRD2). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Detects a ~51 kDa band. Ideal for neuroscience and GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	DRD2 Rabbit pAb
Remarks/Alias	D2R; D2DR; DRD2
Species	Human
Gene ID (Human)	1813
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DRD2 (NP_000786.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14416
Protein Weight	51kDa
Shipping	Ice bag