


Dusp12 Polyclonal Antibody (Rabbit) For Wb And Elisa - Q9Uni6

Item Number: CM0017271



Introduction

DUSP12 Rabbit polyclonal antibody validated for WB and ELISA in human and mouse samples. Detects dual specificity phosphatase 12 (38 kDa) involved in dephosphorylation of tyrosine and serine/threonine residues. For research use.

[Learn More](#)

Parameter	Value
Product Name	DUSP12 Rabbit pAb
Remarks/Alias	YVH1; DUSP1; DUSP12
Species	Human
Gene ID (Human)	11266
Gene ID	11266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-270 of human DUSP12 (NP_009171.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UNI6
Protein Weight	38kDa
Shipping	Ice bag