


## Dusp4 Rabbit Polyclonal Antibody For Wb, Elisa - Q13115

Item Number: CM0026145



### Introduction

DUSP4 rabbit polyclonal antibody detects human DUSP4/MKP-2 in WB and ELISA. Targets the dual-specificity phosphatase that dephosphorylates ERK1/2, with a predicted molecular weight of ~43 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | DUSP4 Rabbit pAb                                                                                                 |
| Remarks/Alias    | TYP, HVH2, MKP2, MKP-2, DUSP4                                                                                    |
| Species          | Human                                                                                                            |
| GeneID (human)   | 1846                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-168 of human DUSP4 (NP_001385.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | Q13115                                                                                                           |
| Protein Weight   | 43kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |