

Dut Rabbit Mab For Wb, Elisa - P33316

Item Number: CM0006605



Introduction

High-quality rabbit monoclonal antibody targeting human DUT (dUTPase), validated in WB and ELISA. Recombinant C-terminal immunogen; 27 kDa; cross-reacts with human. Ideal for nucleotide metabolism and DNA replication research.

[Learn More](#)

Parameter	Value
Product Name	DUT Rabbit mAb
Remarks/Alias	BMFDMS; dUTPase; DUT
Species	Human
Gene ID (Human)	1854
Gene ID	1854
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 160-252 of human DUT (NP_001020419.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P33316
Protein Weight	27kDa
Shipping	Ice bag