


Dvl2 Rabbit Polyclonal Antibody For If/Icc, Elisa - O14641

Item Number: CM0027382



Introduction

DVL2 Rabbit Polyclonal Antibody validated for IF/ICC and ELISA, with cross-reactivity to Mouse. Targets human DVL2 (UniProt O14641), a key mediator of Wnt signaling. Ideal for developmental and signaling research.

[Learn More](#)

Parameter	Value
Product Name	DVL2 Rabbit pAb
Remarks/Alias	DVL2
Species	Human
GeneID (Human)	1856
GeneID	1856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 600-700 of human DVL2 (NP_004413.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	O14641
Protein Weight	79kDa
Shipping	Ice bag