


Dync1Li2 Rabbit Polyclonal Antibody For Wb - 043237

Item Number: CM0021933



Introduction

Rabbit polyclonal antibody to dynein light intermediate chain 2 (DYNC1LI2), involved in microtubule-based vesicle trafficking; validated for WB, IHC-P, ELISA in human, mouse, and rat.

[Learn More](#)

Parameter	Specification
Product Name	DYNC1LI2 Rabbit Polyclonal Antibody
Remarks/Alias	LIC2; DNCL12; DYNC1LI2
Species	Human
Gene ID (Human)	1783
Gene ID	1783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 328-492 of human DYNC1LI2 (NP_006132.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43237
Protein Weight	54kDa
Shipping	Ice bag