

Dysferlin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - 075923

Item Number: CM0003779



Introduction

Dysferlin Rabbit Monoclonal Antibody validated for Western blotting, IHC-P, and ELISA. Specific to human Dysferlin (237 kDa), a calcium sensor in membrane repair. Ideal for neuromuscular disease and muscle biology research applications.

[Learn More](#)

Parameter	Value
Product Name	Dysferlin Rabbit mAb
Remarks/Alias	MMD1; FER1L1; LGMD2B; LGMDR2; Dysferlin (Romeo)
Species	Human
GeneID (Human)	8291
GeneID	8291
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Dysferlin (Romeo) (075923).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	075923
Protein Weight	237kDa
Shipping	Ice bag