

E2F6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa, Chip - O75461

Item Number: CM0027995

Introduction

Rabbit polyclonal antibody against human E2F6 transcription factor. Suitable for WB, IF/ICC, IP, ELISA, and ChIP applications. Cross-reacts with human and mouse. Useful for cell cycle and gene regulation research.

[Learn More](#)

Parameter	Value
Product Name	E2F6 Rabbit pAb
Remarks/Alias	E2F-6; E2F6
Species	Human
GeneID	1876
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-281 of human E2F6 (NP_937987.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse
SWISS	O75461
Protein Weight	32kDa
Shipping	Ice bag