


## Ebf1 Rabbit Pab For Wb, Elisa - Q9Uh73

Item Number: CM0016207



### Introduction

High-quality rabbit polyclonal antibody against human EBF1, validated in WB and ELISA. Ideal for B-cell research and IVD assay development, targeting the pioneer transcription factor COE1.

[Learn More](#)

| Parameter        | Value                                                                                                |
|------------------|------------------------------------------------------------------------------------------------------|
| Product Name     | EBF1 Rabbit pAb                                                                                      |
| Remarks/Alias    | EBF; COE1; OLF1; O/E-1; EBF1                                                                         |
| Species          | Human                                                                                                |
| GeneID           | 1879                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 500-591 of human EBF1 (NP_076870.1) |
| Source           | Rabbit                                                                                               |
| Category         | Polyclonal Antibodies                                                                                |
| Application      | WB, ELISA                                                                                            |
| Cross Reactivity | Human                                                                                                |
| SWISS            | Q9UH73                                                                                               |
| Protein Weight   | 64kDa                                                                                                |
| Shipping         | Ice bag                                                                                              |