

Ece1 Rabbit Pab - P42892

Item Number: CM0019453



Introduction

Rabbit polyclonal antibody targeting human ECE1 (endothelin-converting enzyme 1). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects an ~87 kDa band. Ideal for endothelin-1 processing research.

[Learn More](#)

Parameter	Value
Product Name	ECE1 Rabbit pAb
Remarks/Alias	ECE; ECE1
Species	Human
GeneID (Human)	1889
GeneID	1889
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 670-770 of human ECE1 (NP_001388.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42892
Protein Weight	87kDa
Shipping	Ice bag