


Ech1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13011

Item Number: CM0013466



Introduction

ECH1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, rat ECH1; targets delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase in fatty acid beta-oxidation; for mitochondrial and peroxisomal research.

[Learn More](#)

Parameter	Value
Product Name	ECH1 Rabbit pAb
Remarks/Alias	HPXEL; ECH1
Species	Human
GeneID (Human)	1891
GeneID	1891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-328 of human ECH1 (NP_001389.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13011
Protein Weight	36kDa
Shipping	Ice bag