

Eda Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q92838

Item Number: CM0026848



Introduction

High-quality rabbit polyclonal antibody targeting human EDA (Ectodysplasin-A). Suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reactivity with mouse and rat. Research use in developmental biology and ectodermal dysplasia studies.

[Learn More](#)

Parameter	Description
Product Name	EDA Rabbit pAb
Remarks/Alias	ED1; HED; EDA1; EDA2; HED1; ODT1; XHED; ECTD1; XLHED; ED1-A1; ED1-A2; EDA-A1; EDA-A2; TNLG7C; STHAGX1; EDA
Species	Human
Gene ID (Human)	1896
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 62-182 of human EDA (NP_001390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92838
Protein Weight	41kDa
Shipping	Ice bag