

Eed Rabbit Polyclonal Antibody For Wb, If, Ip, Elisa - 075530

Item Number: CM0027655



Introduction

Rabbit polyclonal antibody against human EED (Polycomb protein EED, HEED, WAIT-1). Suitable for WB, IF/ICC, IP, and ELISA; cross-reacts with human, mouse, and rat. Ideal for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	EED Rabbit pAb
Remarks/Alias	HEED; COG15; WAIT1; EED
Species	Human
Gene ID (Human)	8726
Gene ID	8726
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-400 of human EED (NP_694536.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075530
Protein Weight	50kDa
Shipping	Ice bag