


Eef2 Rabbit Mab - P13639

Item Number: CM0006228



Introduction

Recombinant monoclonal antibody targeting human elongation factor 2 (eEF2). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for translation elongation and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	eEF2 Rabbit mAb
Remarks/Alias	EF2; EF-2; EEF-2; SCA26; eEF2
Species	Human
GeneID (Human)	1938
GeneID	1938
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 759-858 of human eEF2 (P13639).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13639
Protein Weight	95kDa
Shipping	Ice bag