


## Eefsec Rabbit Pab - P57772

Item Number: CM0019514



### Introduction

Rabbit polyclonal antibody against human EEFSEC (SELB), validated for Western blot (WB) and ELISA. Recognizes the selenocysteine-specific elongation factor involved in selenoprotein translation. Useful for studying selenoprotein biosynthesis.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | EEFSEC Rabbit pAb                                                                                                   |
| Remarks/Alias    | SELB, EFSEC, EEFSEC                                                                                                 |
| Species          | Human                                                                                                               |
| GeneID (human)   | 60678                                                                                                               |
| GeneID           | 60678                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 347-596 of human EEFSEC (NP_068756.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | P57772                                                                                                              |
| Protein Weight   | 65kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |