

Egf Receptor (E746-A750Del Specific) Rabbit Polyclonal Antibody For Wb, Elisa - P00533

Item Number: CM0018931



Introduction

Rabbit polyclonal antibody detects EGF receptor with E746-A750 deletion. Validated for WB and ELISA, reacts with human samples. Suitable for cancer research and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	EGF Receptor (E746-A750del Specific) Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant protein of Human EGF Receptor (E746-A750del Specific)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag