


Egfr Rabbit Mab For Wb, Ihc-P, Elisa - P00533

Item Number: CM0007424



Introduction

High-affinity anti-EGFR rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Supports research on EGF receptor signaling, cancer biology, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	EGFR Rabbit mAb
Remarks/Alias	EGFR; ERBB; ERBB1; HER1; NISBD2; PIG61; mENA; epidermal growth factor receptor
Species	Human
Gene ID	1956
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag