


Egfr Rabbit Pab For Wb, Ihc-P, If/Icc, Ip, Elisa - P00533

Item Number: CM0016303



Introduction

EGFR Rabbit pAb (polyclonal) validated for WB, IHC-P, IF/ICC, IP, ELISA applications. Cross-reacts with human, mouse, and rat. Recombinant immunogen spanning amino acids 1021-1210. Suitable for cancer and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	EGFR Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; EGFR
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1021-1210 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag