


Egfr Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P00533

Item Number: CM0020151



Introduction

Rabbit polyclonal antibody targeting human EGFR (ErbB-1, HER1). Validated for WB, IP, ELISA; cross-reacts with human, mouse, rat. Ideal for studying receptor tyrosine kinase signaling and oncology. 134 kDa.

[Learn More](#)

Parameter	Value
Product Name	EGFR Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; EGFR
Species	Human
GeneID	1956
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1100-1200 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag