

Ehd4 Rabbit Pab For Wb And Elisa - Q9H223

Item Number: CM0014633

Introduction

Rabbit polyclonal antibody targeting EHD4/PAST4, validated for Western blot and ELISA. Detects human and mouse EHD4 (~61kDa). EHD4 functions in early endosomal transport and angiogenesis, relevant for cancer and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	EHD4 Rabbit pAb
Remarks/Alias	PAST4; EHD4
Species	Human
Gene ID	30844
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 362-541 of human EHD4 (NP_644670.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H223
Protein Weight	61kDa
Shipping	Ice bag