


Eif3F Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, If/Icc, Ip, Elisa - O00303

Item Number: CM0028945



Introduction

Rabbit polyclonal antibody against human EIF3F (eIF3-p47, MRT67, EIF3S5). Validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat. Suitable for translation initiation and Notch signaling research. ~38 kDa.

[Learn More](#)

Parameter	Value
Product Name	EIF3F Rabbit pAb
Remarks/Alias	MRT67; EIF3S5; eIF3-p47; EIF3F
Species	Human
GeneID (human)	8665
GeneID	8665
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 88-357 of human EIF3F (NP_003745.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00303
Protein Weight	38kDa
Shipping	Ice bag