


Eif3H Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Elisa - O15372

Item Number: CM0014989



Introduction

Rabbit polyclonal antibody against EIF3H (eIF3-p40), validated for WB, IP, IF/ICC, ELISA in human, mouse, rat. Ideal for studying translation initiation and protein synthesis. - O15372

[Learn More](#)

Parameter	Value
Product Name	EIF3H Rabbit pAb
Remarks/Alias	EIF3S3; eIF3-p40; eIF3-gamma; EIF3H
Species	Human
GeneID (Human)	8667
GeneID	8667
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-340 of human EIF3H (NP_003747.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15372
Protein Weight	40kDa
Shipping	Ice bag