

Eif4Enif1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nra8

Item Number: CM0017659



Introduction

EIF4ENIF1 Rabbit Polyclonal Antibody (pAb) for WB, IF/ICC, ELISA. Targets eukaryotic translation initiation factor 4E transporter, a key P-body component regulating mRNA stability and translation. Cross-reacts with human, mouse, rat. Ideal for translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	EIF4ENIF1 Rabbit pAb
Remarks/Alias	4E-T; Clast4; EIF4ENIF1
Species	Human
GeneID	56478
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human EIF4ENIF1 (NP_001157973.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRA8
Protein Weight	108kDa
Shipping	Ice bag