

Eif4G2/P97 Rabbit Mab For Wb, Ihc-P, Ip, And Elisa - P78344

Item Number: CM0006275



Introduction

EIF4G2/p97 Rabbit Monoclonal Antibody for western blot, immunohistochemistry, immunoprecipitation, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for translation initiation and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	EIF4G2/p97 Rabbit mAb
Remarks/Alias	P97; AAG1; DAP5; NAT1; EIF4G2/p97
Species	Human
GeneID (human)	1982
GeneID	1982
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-490 of human EIF4G2/p97 (NP_001409.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78344
Protein Weight	102kDa
Shipping	Ice bag