

Eif5A Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P63241

Item Number: CM0023160



Introduction

Polyclonal antibody targeting human eIF5A, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse eIF5A, involved in translation elongation and ribosome rescue. Ideal for research on protein synthesis, apoptosis, and viral replication.

[Learn More](#)

Parameter	Value
Product Name	eIF5A Rabbit pAb
Remarks/Alias	FABAS; EIF-5A; EIF5A1; eIF-4D; eIF5A1; eIF5A
Species	Human
Gene ID (Human)	1984
Gene ID	1984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human eIF5A (NP_001137234.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P63241
Protein Weight	17kDa
Shipping	Ice bag