


Elf6 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P56537

Item Number: CM0022452



Introduction

Rabbit polyclonal antibody against human EIF6, validated for WB, IF/ICC, ELISA. Detects endogenous EIF6 (27kDa) with cross-reactivity to mouse and rat. Suitable for research on translation initiation and ribosome biogenesis.

[Learn More](#)

Parameter	Value
Product Name	EIF6 Rabbit pAb
Remarks/Alias	CAB; EIF3A; eIF-6; p27BBP; ITGB4BP; b(2)gcn; p27(BBP); EIF6
Species	Human
Gene ID	3692
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human EIF6 (NP_852133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56537
Protein Weight	27kDa
Shipping	Ice bag