

Elac2 Rabbit Pab - Q9Bq52

Item Number: CM0028690



Introduction

ELAC2 rabbit polyclonal antibody validated for Western blot (WB), immunohistochemistry (IHC-P), and ELISA; detects human, mouse, and rat ELAC2, a mitochondrial tRNA processing endonuclease involved in RNA maturation and ribosome assembly.

[Learn More](#)

Parameter	Value
Product Name	ELAC2 Rabbit pAb
Remarks/Alias	ELC2; HPC2; COXPD17; ELAC2
Species	Human
GeneID (Human)	60528
GeneID	60528
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQ52
Protein Weight	92kDa
Shipping	Ice bag