

Elapor1 Rabbit Pab - Q6Uxg2

Item Number: CM0020119

Introduction

ELAPOR1 Rabbit pAb is a polyclonal antibody for WB and ELISA, detecting human, mouse, and rat ELAPOR1 (111kDa). This endosome/lysosome-associated apoptosis and autophagy regulator is ideal for cancer and cell death research.

[Learn More](#)

Parameter	Value
Product Name	ELAPOR1 Rabbit pAb
Remarks/Alias	EIG121; KIAA1324; ELAPOR1
Species	Human
GeneID (human)	57535
GeneID	57535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human ELAPOR1 (NP_065826.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UXG2
Protein Weight	111kDa
Shipping	Ice bag