

Elmod3 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Fg2

Item Number: CM0026122

Introduction

ELMOD3 rabbit polyclonal antibody validated for WB, ELISA, reacts with human, mouse, rat. Target is a GTPase-activating protein for ARL2, implicated in hearing loss. Ideal for IVD and research.

[Learn More](#)

Parameter	Value
Product Name	ELMOD3 Rabbit pAb
Remarks/Alias	LST3; RBED1; RBM29; DFNA81; DFN888; ELMOD3
Species	Human
GeneID (Human)	84173
GeneID	84173
Immunogen	Recombinant fusion protein corresponding to amino acids 290-381 of human ELMOD3 (NP_001128493.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FG2
Protein Weight	43kDa
Shipping	Ice bag