

Elovl6 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H5J4

Item Number: CM0020793

Introduction

High-quality rabbit polyclonal antibody targeting human ELOVL6, suitable for Western blotting and ELISA. Detects ELOVL6 in human, mouse, and rat samples. Ideal for studies on very long chain fatty acid elongation and lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	ELOVL6 Rabbit pAb
Remarks/Alias	FAE; LCE; FACE; hELO2; ELOVL6
Species	Human
Gene ID	79071
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human ELOVL6 (NP_076995.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H5J4
Protein Weight	31kDa
Shipping	Ice bag