

Eno1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P06733

Item Number: CM0015055



Introduction

Rabbit polyclonal antibody targeting human α -enolase (ENO1) for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse. Ideal for studies on glycolysis, transcriptional repression, and plasminogen activation.

[Learn More](#)

Parameter	Value
Product Name	ENO1 Rabbit pAb
Remarks/Alias	NNE; PPH; MPB1; ENO1L1; ENO1-IT1; HEL-S-17; ENO1
Species	Human
GeneID (Human)	2023
GeneID	2023
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-434 of human ENO1 (NP_001419.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P06733
Protein Weight	47kDa
Shipping	Ice bag