

Eno2 Rabbit Mab For Wb - P09104

Item Number: CM0002804



Introduction

Rabbit monoclonal antibody targeting human ENO2 (gamma-enolase, NSE), validated for WB, IF/ICC, IF-F, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience research and neuron-specific enolase detection.

[Learn More](#)

Parameter	Value
Product Name	ENO2 Rabbit mAb
Remarks/Alias	NSE; HEL-S-279; ENO2
Species	Human
GeneID (Human)	2026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 47-100 of human ENO2 (NP_001966.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09104
Protein Weight	47 kDa
Shipping	Ice bag