

Eno2 Rabbit Polyclonal Antibody For Wb And Elisa - P09104

Item Number: CM0026027



Introduction

ENO2 Rabbit Polyclonal Antibody targeting human gamma-enolase (NSE) for Western blot and ELISA. Validated cross-reactivity with human, mouse, and rat. Ideal for studying glycolysis, neuroprotection, and cancer biomarkers.

[Learn More](#)

Specification	Details
Product Name	ENO2 Rabbit pAb
Remarks/Alias	NSE; HEL-S-279; ENO2
Species	Human
Gene ID (Human)	2026
Gene ID	2026
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human ENO2 (NP_001966.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09104
Protein Weight	47kDa
Shipping	Ice bag