

Enos Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P29474

Item Number: CM0013675



Introduction

High-quality rabbit polyclonal antibody against human endothelial nitric oxide synthase (eNOS/NOS3). Validated for WB, IF-P, ELISA. Detects 133 kDa band; cross-reacts with mouse and rat. Ideal for vascular and NO signaling research.

[Learn More](#)

Parameter	Value
Product Name	eNOS Rabbit pAb
Remarks/Alias	eNOS; ECNOS
Species	Human
GeneID (Human)	4846
GeneID	4846
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1000-1203 of human eNOS (NP_000594.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29474
Protein Weight	133kDa
Shipping	Ice bag