

## Enpep Rabbit Pab For Wb, Elisa - Q07075

Item Number: CM0012858



### Introduction

ENPEP Rabbit pAb (APA/CD249/gp160) is a polyclonal antibody detecting human and mouse glutamyl aminopeptidase in WB and ELISA. Ideal for hypertension and renal research.

[Learn More](#)

| Parameter        | Value                    |
|------------------|--------------------------|
| Product Name     | ENPEP Rabbit pAb         |
| Remarks/Alias    | APA; CD249; gp160; ENPEP |
| Species          | Human                    |
| GeneID (human)   | 2028                     |
| GeneID           | 2028                     |
| Immunogen        | Synthetic peptide        |
| Source           | Rabbit                   |
| Category         | Polyclonal Antibodies    |
| Application      | WB, ELISA                |
| Cross Reactivity | Human, Mouse             |
| SWISS            | Q07075                   |
| Protein Weight   | 109kDa                   |
| Shipping         | Ice bag                  |