

Enpep Rabbit Pab For Wb, Ihc-P, Elisa - Q07075

Item Number: CM0026946



Introduction

Accurate detection of glutamyl aminopeptidase across human, mouse, and rat samples. This rabbit pAb supports Western blot, IHC-P, and ELISA in hypertension and renal biology.

[Learn More](#)

Parameter	Value
Product Name	ENPEP Rabbit pAb
Remarks/Alias	APA; CD249; gp160; ENPEP
Species	Human
GeneID (Human)	2028
GeneID	2028
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 708-957 of human ENPEP (NP_001968.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07075
Protein Weight	109kDa
Shipping	Ice bag