


Enpp2 Rabbit Pab For Wb, Elisa - Q13822

Item Number: CM0020502



Introduction

ENPP2 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, targeting human and mouse ENPP2/ATX lysophospholipase D (Q13822) for lysophosphatidic acid pathway research.

[Learn More](#)

Parameter	Value
Product Name	ENPP2 Rabbit pAb
Remarks/Alias	ATX; NPP2; ATX-X; PDNP2; LysoPLD; AUTOTAXIN; PD-IALPHA; ENPP2
Species	Human
GeneID (Human)	5168
GeneID	5168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 564-863 of human ENPP2 (NP_001035181.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13822
Protein Weight	99kDa
Shipping	Ice bag